IN NURSING HOMES

PDA SOLUTIONS




HEB'DOSE DUO DISPENSER

THE ORIGINAL SOLUTION


A solution which has been on the market for over 10 years and has proved its worth to the professional user, in nursing homes.

HEB'DOSE © DUO PDA

- A3 format weekly dispenser composed of 28 double compartments, detachable by the day and by dispensing time.
- An exceptional compartment to easily manage special doses, additional treatments or crushable / non crushable medications.
- A two-step opening to comply with dosing recommendations.
- Perfectly integrated in the nursing home environment or at home, thanks to an adapted distribution system.




PRODUCTS' HIGHLIGHTS COMMON TO THE ENTIRE RANGE



HEB'DOSE PDA ADAPTATED SOLUTIONS FOR NURSING HOMES

SECUR'DOSE[©] SOFTWARE

TO TRACE AND SECURE THE MEDICATION CIRCUIT FROM PREPARATION TO PATIENT ADMINISTRATION



HEB'DOSE RANGE PILL CONTAINERS

ERGONOMIC MATERIAL ACCESSIBLE AND UNDERSTANDABLE BY ALL

The single use dispensers are eco-designed, guaranteeing compliance with hygiene and safety standards. In APET food contact, they are composed of virgin material, in contact with the drug, and regenerated material. Produced in France in a cleanroom environment, they are recyclable.

THE OPENING



THEY HAVE BEEN DESIGNED TO FACILITATE : THE TRANSMISSION

TAKING THE MEDICATION



md se • Tél. +33 (0)1 45 10 77 90 • Fax +33 (0)1 45 10 15 95 • www.mdose.fr

HEB'DOSE RANGE PILL CONTAINERS

GUIDED BY THE SÉCUR'DOSE SOFTWARE

To adapt to the caregiver's expectations and constraints.


HEB'DOSE[©] SOLO PDA

- A4 format weekly dispenser composed of 28 single compartments, detachable by the day and by dispensing time.
- A different capacity for dispensing times, adapted to the average number of medications usually prescribed.
- Ideal for establishments managing patients with a light treatment of for counter patients.
- Economical solution.

Development Biological accession Biologica


PDA HEB'DOSE © MULTI

- A3 format weekly dispenser composed of 28 multicompartments per dispensing time, detachable by the day and by dispensing time.
- A different number of compartments by dispensing time to optimize the pill container's capacity and avoid waste.
- A sealing operation carried out by a single label that includes the essential information.
- Ideal for establishments wishing to separate medications from others.



HEB'DOSE © XL PDA

- A3 format weekly dispenser composed of 28 XL compartments, detachable by the day and by dispensing time.
- 4 XL compartments of different capacities, adapted to the average number of medications usually prescribed.
- A one-step sealing operation.
- Perfect integration into our distribution solutions in nursing homes.
- For establishments wishing to keep primary packaging.



DISTRIBUTION SOLUTIONS

IN NURSING HOMES


To meet safety and traceability requirements, to limit administration errors and facilitate the follow-up.

«PLUMIERS» FOR NURSING HOME DISTRIBUTION

SIMPLIFY AND SECURE THE CAREGIVERS' WORK

- A daily «plumier», one patient, 6 identification criteria.
- An additional nominative storage for sachets, primary unit packs or treatment additions.
- Perfectly smooth for easy cleaning.

- Dispensers positioned on a 600x400 standard tray 8 dispensers per tray.
- Perfect integration into our Caméléon distribution trolleys.



«Plumier» dispenser for Heb'Dose Duo Mdose «Plumier» with a lid for Heb'Dose Multi, XL and Solo

BINDER SOLUTION CLASSEUR

AN OPTIMIZED SCENARIO


- Saves time when preparing the trolley.
- Binders can contain up to 16 weekly PDA depending on the dispenser.
- Pre-cuts which facilitate the picking of each dispensing time and avoid an excessive handling of the blister.
- Off-dispenser and in 600x400 baskets storage solution.
- Perfect integration into our Caméléon trolleys.





DISTRIBUTION ESSENTIALS

IN NURSING HOMES





Julitet 2023 N° d'enregistrement en application de l'article L. 541-10 du code de l'Environnement : FR320613_01SSIS / FR320613_03QMQX